



Dentek Temperin Max Repair Kit

Uses

Temporary treatment to replace a lost filling, cap, crown or inlay for up to 48 hours until a dentist is seen.

Dosage

Follow step-by-step instructions on the pack. A dentist must be visited within 48 hours of use.

Format and Ingredients

Dental repair kit containing:
Zinc Oxide, Calcium Sulfate, Paraffinum Liquidum, Potassium Alum, Aluminium Sulfate, Aluminium Phosphate, Menthol, Eugenol, Iron Oxide Yellow,

Category:
Mouth Care, Analgesics

Manufacturer:
PRESTIGE CONSUMER HEALTH

Pack size:

RRP: £4.99

Legal Status: MD